

High Sensitivity Fluorescence Detector for Supercritical Fluid Chromatography

Satoe Iijima¹, Akitaka Terada¹, DJ Tognarelli², John Burchell², Takeshi Kanomata¹, Masao Bounoshita¹, Yasuyo Sato¹, Miki Kuwajima¹

¹JASCO Corporation, 2967-5 Ishikawa-machi, Hachioji, Tokyo 192-8537

²JASCO Incorporated, 28600 Mary's Court, Easton, MD 21601

E-mail: satoe.iijima@jasco.co.jp

Introduction

Fluorescence detection has been widely used as one of a practical detection method in HPLC analysis. It enables the selective detection of fluorescent substances, with sensitivity that can be up to 1000 times greater than with UV detection.

Fluorescence detection can also be applied to non-fluorescent substances with the use of derivatization, with many derivatization agents becoming commercially available in recent years. Therefore, this technique provides the advantages of highly sensitivity and selectivity to many different classes of compound, and broadens the range of applications in many industries. However, fluorescence detection has not been used in SFC analysis due to the difficulty of developing a suitable high pressure flow cell.

We developed a fluorescence detector with a flow cell specifically constructed for operations at the high pressures required for SFC. This fluorescence detector is included in our newer SFC system (SFC-4000 series). In this presentation, we will introduce an application for the analysis of several compounds using SFC separation with fluorescence detection.

Experimental

Apparatus

Figure 1 shows the JASCO SFC system equipped with UV/Visible (UV) and fluorescence (FL) detectors used in this experiment. Figure 2 also shows the schematic diagram of this system. This system enables method scouting analysis to be performed on several mobile phases and columns.

Figure 3 shows the structures and differences between a conventional and novel fluorescence detector flow cell. We developed a unique quartz flow cell design and structure for use with high pressure (up to 20 MPa). The outer surfaces of this flow cell are coated with aluminum, and provides high collection efficiency of fluorescence by functioning as a reflector.








Fig.3 Structures and differences between the conventional and novel fluorescence flow cell.

Results and Discussion


SFC/UV vs. SFC/FL




SFC/FL Applications



Vitamin E (α, β, γ, δ-Tocopherol)



Quinolone Antibiotic



Conclusion

- We developed a fluorescence detector with a flow cell specifically constructed for operation at the high pressures required for SFC
- SFC/FL provides highly sensitive detection (from fg to pg) for fluorescent substances, compared to UV/Visible detector.
- SFC using fluorescence detection will open the door to new applications fields.

References

- 1.V. R. Meyer, "Practical High-Performance Liquid Chromatography", 4th ed., John Wiley & Sons, 2004, 92 and 279-282.
- 2.M. Bounoshita, T. Kanomata, CHROMATOGRAPHY, 2011, 32, 23-32.
- 3.O. A. Adegoke, Afr.J. Pure Appl.Chem., 2012, 6 (14), 129-140.
- 4.T. Yarita, A. Nomura, K. Abe, Y. Takeshita, J. Chromatogr. A, 1994, 679, 329-334.