

Analyzing the interaction of human serum albumin and 3,5-diiodosalicylic acid probed by circular dichroism

Introduction

Human serum albumin (HSA) is the most abundant protein in blood plasma. HSA binds with pharmaceutical compounds and other in-vivo substances and plays an important role in the transport of these substances to target organs. Several studies have reported on the binding affinity of HSA and its interactions with a variety of compounds.¹⁻³

This application notes demonstrates the use of the J-1500 CD spectrometer and ATS-530 automatic titrator to monitor the interaction of HSA titrated with 3,5-diiodosalicylic acid. While 3,5-diiodosalicylic acid is achiral, its interaction with chiral HSA induces circular dichroism.

Keywords

J-1500, circular dichroism, HSA, pharmaceuticals, ATS-530 automatic titrator, biochemistry, pharmaceutical

Experimental

Measurement conditions	
Data acquisition interval	2 sec
Scan speed	100 nm/min
Spectral bandwidth	1 nm
Data pitch	0.1 nm
Path length	10 mm
Accumulations	2



JASCO J-1500 CD spectrometer
View product information at www.jascoinc.com

50 mL aliquots of 0.025 mM of 3,5-diiodosalicylic acid was added to 2 mL of 0.0228 mM of HSA in 100 mM acetate buffer (pH 6.3) 20 times. The CD spectrum was scanned from 360 to 260 nm and also monitored at 320 nm.

Results

Chirality can be induced in an achiral substance interacting with a chiral substance. This interaction will exhibit circular dichroism. The resulting CD spectrum from the interaction of achiral 3,5-diiodosalicylic acid with chiral HSA shows a positive peak at 320 nm in Figure 1. While HSA does not show a CD signal at 320 nm, the induced CD from the interaction of HSA and 3,5-diiodosalicylic acid depicts this signal at 320 nm which increases with the increasing additions of 3,5-diiodosalicylic acid and is shown in Figure 2.




Figure 1. 3-D graphical representation of the CD spectrum of HSA titrated with 3,5-diiodosalicylic acid.




Figure 2. CD spectra of HSA titrated with 3,5-diiodosalicylic acid. The arrows represent the increasing volume of 3,5-diiodosalicylic acid.

The increase in the induced CD signal at 320 nm was plotted in Figure 3 as a function of the 3,5-diiodosalicylic acid concentration. Figure 3 is also known as a Hill plot which describes the cooperativity of the interaction between HSA and 3,5-diiodosalicylic acid. A dissociation constant, which is the concentration at which 50% of a substance is bound, can also be estimated from a Hill plot. Using Figure 3, the dissociation constant of 3,5-diiodosalicylic acid binding to HSA was determined to be 0.023 mM. The Hill coefficient is approximately 3.1 which indicates a positive cooperative reaction.




Figure 3. CD signal monitored at 320 nm for the titration of 3,5-diiodosalicylic acid to HSA.

Conclusion

This application note demonstrates the use of the ATS-530 automatic titrator and J-1500 CD spectrometer to obtain CD spectra of the titration of 3,5-diiodosalicylic to HSA in order to determine a dissociation constant.

References

1. R. Brodersen, J. Biol. Chem. (1977) 252, 14, 5067-5072.
2. U. Kragh-Hansen, Biochem. J. (1983) 209, 135-142.
3. Sinha, S. S., Mitra, R. K., and S. K. Pal, J. Phys. Chem.B. (2008) 112, 4884-4891.